



## HOW TO USE

Designed for  
horses and the  
people who care  
for them

# How to use Aservo<sup>®</sup> EquiHaler<sup>®</sup>



**Breath Indicator** - the timing of actuation should match the beginning of inspiration, seen as an inward deflection of the BreathView<sup>®</sup> membrane

**Air Inlet** - source of air

**Handle** - prior to first use, the Aservo<sup>®</sup> EquiHaler<sup>®</sup> must be activated and primed. Activation consists of piercing the cartridge and priming consists of 3 actuations until the Soft Mist<sup>™</sup> appears

**Piercing Element with Fill Indicator** - the fill indicator provides guidance for the user on the volume of medication remaining in the inhaler



## 10 Days total treatment duration

<b>DAY 1-5</b>	<b>8 ACTUATIONS (2,744 µg) TWICE DAILY</b>
<b>DAY 6-10</b>	<b>12 ACTUATIONS (4,116 µg) ONCE DAILY</b>

**Nostril Adaptor** - is inserted into the left nostril

**Lever** - Actuation is a 2 step process that consists of a full pull on the lever (to load the spring), and a gentle pull on the lever (to release the mist)